

NI LabView Installation Instructions

Chrissy Willow Rothgeb | Cise-IT@jmu.edu | updated 2024-04-05

Download Instructions

Go to [https://support.cise.jmu.edu/installs/NI/Web Installers/](https://support.cise.jmu.edu/installs/NI/Web%20Installers/) and log in with your JMU eID.

support.cise.jmu.edu - /installs/NI/Web Installers/

[\[To Parent Directory\]](#)

Monday, January 8, 2024	1:45 PM	1069582004	LabVIEW2023Q3Full.dmg
Tuesday, July 12, 2022	8:34 AM	6769400	ni-software-platform-bundle_21.5_suite_online.exe
Monday, January 8, 2024	1:33 PM	4897464	ni-software-platform-bundle_23.5_suite_online.exe
Tuesday, February 6, 2024	10:02 AM	<dir>	old

Windows users will need to download and run the Web Installer: [ni-software-platform-bundle 21.5 suite online.exe](#). This small application will allow you to install any application from the LabView suite of programs via the internet.

Note: NI LabView 2023 is not compatible with the MyDAQ and it is advised you do not use this version if you plan to use MyDAQ hardware.


Mac users will need to download and run the DMG installer: [LabVIEW2023Q3Full.dmg](#). This includes individual installs for LabView 2021, and LabView 2021 Runtime Engine.

Windows Installation Process

NOTE: If you already have LabView Installed, please uninstall all prior products before proceeding!

- <https://www.ni.com/.../uninstall-national-instruments-software-or-drivers-using-ni-pack.html>

Open the [web installer](#) you downloaded. The Software Platform Bundle will install the NI Package Manager before continuing with LabView. You may choose to disable Fast Windows Startup if you wish. By default, the installer will select LabVIEW, NI-DAQmx, NI-488.2, NI-VISA, PXI Platform Services, and NI-Sync. At a minimum, you will also want to **add DAQExpress, NI-Sync, and NI-ELVISmx.**



You may select as many products as you like, but please be aware, installing *all* packages may take over 80GB of space. **The minimum install will take roughly 14GB.** You can always run the web installer later to add any additional products.

Accept the license agreements and proceed with the install. This may take a half hour or more depending upon your product selection.

After installation has completed, you will be asked to log in with an NI Account. If you do not already have an account, click [Create a user account](#) from the installer. Once you have logged in, choose to **Enter a serial number**, and **type the 9-digit number provided by your professor** into each of the fields.


The screenshot shows a window titled "Activate Software" with a close button (X) in the top right corner. The user is logged in as "rothgecw@jmu.edu" with a "Log Out" link. The main heading is "7 licenses selected". Below this, a message states: "A serial number is a unique identifier belonging to each product that allows online product activation." with a "Learn more" link. There is a text input field labeled "Enter a serial number" with a dropdown arrow and an "Activate" button. Below the input field is a list of products, each with a text input field for a serial number: "LabVIEW 2021 Application Builder", "LabVIEW 2021 Base Development System", "LabVIEW 2021 Debug Deployment System", "LabVIEW 2021 Full Development System", and "LabVIEW 2021 LabVIEW Database Connectivity Toolkit". At the bottom, there is a "Privacy Statement" link and two buttons: "Extend Trial (4)" and "Cancel".

Not all products are covered under student licensing. Click Finish and ignore any remaining activations.

The screenshot shows a window titled "Activate Software" with a close button (X) in the top right corner. The main heading is "Congratulations!". Below this, a message states: "Activation was a success. See below for details." There is a progress indicator showing "4 ACTIVATED" (underlined) and "3 REMAINING". Below the progress indicator is a list of activated products, each with a checkmark and a duration: "223 days - LabVIEW 2021 Application Builder", "223 days - LabVIEW 2021 LabVIEW Database Connectivity Toolkit", "223 days - LabVIEW 2021 Professional Development System", and "223 days - LabVIEW 2021 LabVIEW Report Generation Toolkit for Microsoft Office". At the bottom, there is a "Product Registration" link and two buttons: "Activate Remaining" and "Finish".


Mac Installation Process

The file you downloaded earlier ([LabVIEW2023Q3Full.dmg](#)) may be opened with Finder and contains individual installs for LabVIEW 2021, and LabVIEW 2021 Runtime Engine (under Extras):



NOTE: Driver support on Mac is limited to NI-488.2 and NI-VISA, and must be downloaded separately:
- <https://ni.com/en-us/support/downloads/drivers.html>

If you receive a message that the installer “can’t be opened because Apple cannot check it for malicious software,” open Security & Privacy from the System Preferences app and tell it to Open Anyway:



Support

If you have any difficulty installing or licensing LabView, please email Cise-IT@jmu.edu for support.
Thank you!

Chrissy Willow Rothgeb

Cise IT - Labs and Server Support

College of Integrated Science and Engineering | James Madison University

540.568.2629 | ISAT/CS 260 | 701 Carrier Drive, Harrisonburg VA 22807